

Mouse Monoclonal anti-Cyclin D1, Clone AM29

60-0021; 60-0021-7 61-0021; 61-0021-2; 61-0021-5 Isotype Concentration	6 mL; 7 mL predilute Antibody, Ready-To-Use 1 mL; 0.2 mL; 0.5 mL Concentrate Antibody IgG2a See container label
Intended Use	For In Vitro Diagnostic Use. This product is used to qualitatively detect Cyclin D1 in normal and neoplastic formalin fixed, paraffin embedded (FFPE) tissue sections in immunohistochemical (IHC) detection methodology. Interpretation must be made within the context of the patient's clinical history and other diagnostic test by a qualified pathologist.
Description	Anti-Cyclin D1 antibodies are useful in IHC analysis for the diagnosis of mantle cell lymphomas. It is an important marker for the differentiation between mantle cell lymphomas and other small B-cell lymphomas. In addition, anti-Cyclin D1 antibody is a useful tool to identifying overexpression of cyclin D1 in breast and colon carcinomas.
Reagent provided	This antibody is diluted in 10 mM phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN ₃) as antimicrobial agent.
Precautions	For professional users.
1 roodutione	Proper handling of this product as with any product derived from biological sources according to local and applicable regulations.
	Sodium azide (NaN ₃) is a toxic chemical. The concentration in this product is not classified as hazardous, however, the build-ups of NaN ₃ may react with lead and copper plumbing to form highly explosive metal azides. Flush the disposed reagent with large volume of water to prevent azide build-up.
Usage	
Dilution	60-0021; 60-0021-7: Ready-To-Use
	61-0021; 61-0021-2; 61-0021-5: Dilute 1:50-100 before use when using Acu-Stain [™] detection system. Optimum dilution factor may vary depending on the specimen and preparation process and should be determined by each individual investigator.
Staining procedure	Incubate this antibody with tissue section for 30-60 minutes at room temperature. Follow the instructions from the selected detection system.
Positive control tissue	Mantle Cell Lymphoma or Breast Cancer
Epitope retrieval	HIER, Tris-EDTA pH 9
Staining pattern	Nucleus
Storage	Store at 2-8°C.
References	 Ott M, et al. The Journal of Pathology. 1996; 179(3):238-242 Leong AS-Y, et al. Manual of Diagnostic Antibodies for Immunohistology. 2003; pp. 165-166 de Boer CJ, et al. Blood. 1995 Oct 1;86(7):2715-23.
Symbols	REF LOT IVD Image: Second seco

30021 Rev.02

