

## Mouse Monoclonal anti-Melan-A (MART-1), Clone A103

60-0092; 60-0092-7 6 mL; 7 mL Pre-dilute Antibody, Ready-To-Use

61-0092; 61-0092-2 1 mL; 0.2 mL Concentrate Antibody

Isotype lgG1

Concentration See container label

For In Vitro Diagnostic Use. Intended Use

> This product is used to qualitatively detect Melan-A in normal and neoplastic formalin fixed, paraffin embedded (FFPE) tissue sections in immunohistochemical (IHC) detection methodology. Interpretation must be made within the context of the patient's clinical history and other diagnostic

test by a qualified pathologist.

Melan-A, also known as Melanoma Antigen Recognized by T cells 1 (MART-1), is a melanocyte Description

differentiation antigen. The anti-Melan-A body labels melanocytes, steroid producing cells. It is a useful tool for the identification of melanomas. Melan-A is also expressed in angiomyolipomas.

The antibody is applicable as a useful marker for angiomyolipomas.

This antibody is diluted in 10 mM phosphate buffered saline (PBS), pH 7.2 containing 1% bovine Reagent provided

serum albumin (BSA) and 0.09% sodium azide (NaN<sub>3</sub>) as antimicrobial agent.

**Precautions** For professional users.

Proper handling of this product as with any product derived from biological sources according to

local and applicable regulations.

Sodium azide (NaN<sub>3</sub>) is a toxic chemical. The concentration in this product is not classified as hazardous, however, the build-ups of NaN<sub>3</sub> may react with lead and copper plumbing to form highly explosive metal azides. Flush the disposed reagent with large volume of water to prevent

azide build-up.

Usage

**Dilution** 60-0092; 60-0092-7: Ready-To-Use

61-0092; 61-0092-2: Dilute 1:50-100 before use when using Acu-Stain<sup>™</sup> detection system.

Optimum dilution factor may vary depending on the specimen and preparation process and should

be determined by each individual investigator.

Staining procedure Incubate this antibody with tissue section for 30-60 minutes at room temperature. Follow the

instructions from the selected detection system.

Positive control tissue Melanoma

HIER Tris buffer pH 9 **Epitope retrieval** 

Staining pattern Cytoplasm

Storage Store at 2-8°C.

1. Chen Y-T, et al. Proc Natl Acad Sci 1996;93:5915-9. References

2. Jungbluth AA, et al. Am J Surg Pathol 1998;22:595-602.

**Symbols** 

REF Catalog No.

LOT Batch No.

IVD In Vitro Diagnostic Use



30794 Rev.02





