

## Mouse Monoclonal anti-Calponin, Clone CALP

60-0105; 60-0105-7 61-0105; 61-0105-2; 61-0105-5 Isotype Concentration:	6 mL; 7 mL predilute Antibody, Ready-To-Use 1 mL; 0.2 mL; 0.5 mL Concentrate Antibody IgG1 See container label
Intended Use	For In Vitro Diagnostic Use.
	This product is used to qualitatively detect Calponin in normal and neoplastic formalin fixed, paraffin embedded (FFPE) tissue sections in immunohistochemical (IHC) detection methodology. Interpretation must be made within the context of the patient's clinical history and other diagnostic test by a qualified pathologist.
Description	Calponin interacts with actin, tropomyosin, and calmodulin. It is restrict expressed in smooth muscle cells and is involved in the regulation of smooth muscle contraction. Anti-calponin has been found to be useful in differentiating benign sclerosing lesions of the breast from carcinoma. Calponin positivity has also been noted in malignant myoepithelioma and pleomorphic adenoma of salivary gland origin, as well as angiomatoid malignant fibrous histiocytoma.
Reagent provided	This antibody is diluted in 10 mM phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN $_3$ ) as antimicrobial agent.
Precautions	For professional users.
	Proper handling of this product as with any product derived from biological sources according to local and applicable regulations.
	Sodium azide $(NaN_3)$ is a toxic chemical. The concentration in this product is not classified as hazardous, however, the build-ups of $NaN_3$ may react with lead and copper plumbing to form highly explosive metal azides. Flush the disposed reagent with large volume of water to prevent azide build-up.
Usage	
Dilution	60-0105; 60-0105-7: Ready-To-Use
	61-0105; 61-0105-2; 61-0105-5: Dilute 1:50 to 1:100 before use when using Acu-Stain <sup>™</sup> detection system. Optimum dilution factor may vary depending on the specimen and preparation process and should be determined by each individual investigator.
Staining procedure	Incubate this antibody with tissue section for 30-60 minutes at room temperature. Follow the instructions from the selected detection system.
Positive control tissue	Appendix, breast, leiomyoma, uterus
Epitope retrieval	Proteinase K
Staining pattern	Cytoplasm
Storage	Store at 2-8°C.
References	<ol> <li>Lazard D, et al. Proc Natl Acad Sci USA 1993; 90:999.</li> <li>Gimona M, et al. FEBS Lett 1990; 274(1,2):159</li> </ol>
Symbols	REF     LOT     IVD     Image: Catalog No.     In Vitro Diagnostic Use     Temperature Range     Use By

30809 Rev.02



