

Mouse Monoclonal anti-Prostatic Acid Phosphatase (PsAP or PAP), Clone PASE/4LJ

60-0104; 60-0104-7 6 mL; 7 mL predilute Antibody, Ready-To-Use 61-0104; 61-0104-2; 61-0104-5 1 mL; 0.2 mL; 0.5 mL Concentrate Antibody

Isotype IgG1

Concentration See container label

Intended Use For In Vitro Diagnostic Use.

This product is used to qualitatively detect Prostatic Acid Phosphatase (PsAP or PAP) in normal and neoplastic formalin fixed paraffin embedded (FFPE) tissue sections in immunohistochemical detection methodology. Interpretation must be made within the context of the patient's clinical

history and other diagnostic test by a qualified pathologist.

Description The antibody labels glandular epithelium of the prostate and is a useful tool for the identification of

normal prostatic tissue and prostatic carcinoma.

Reagent provided This antibody is diluted in 10 mM Phosphate buffered saline (PBS), pH 7.2 containing 1% bovine

serum albumin (BSA) and 0.09% sodium azide (NaN₃) as antimicrobial agent.

Precautions For professional users.

Proper handling of this product as with any product derived from biological sources according to

local and applicable regulations.

Sodium azide is a toxic chemical. The concentration in this product is not classified as hazardous, however, the build-ups of NaN_3 may react with lead and copper plumbing to form highly explosive metal azides. Flush the disposed reagent with large volume of water to prevent azide build-up.

Usage

Dilution 60-0104; 60-0104-7: Ready-To-Use

61-0104; 61-0104-2; 61-0104-5: Dilute 1:50-100 before use when using Acu-Stain[™] detection system. Optimum dilution factor may vary depending on the specimen and preparation process

and should be determined by each individual investigator.

Staining procedure Incubate this antibody with tissue section for 30-60 minutes at room temperature. Follow the

instructions from the selected detection system.

Positive control tissue Prostate

Epitope retrieval Not Required

Staining pattern Cytoplasm

Storage Store at 2-8°C.

References 1. Haines A, et al. Biochem Soc Trans 1987;15:1179-80.

2. Haines A, et al. Br J Cancer 1989;60:887-92.

3. Nadji M, et al. Am J Clin Path 1980;73:735-9.

Symbols

REF
Catalog No.

LOT Batch No. In Vitro Diagnostic Use

Temperature Range





